Carl Zeiss do Brasil Ltda.

ISO 9001:2000
Reg. No. 311013 QM

Carl Zeiss do Brasil

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Tel: 11

Fax: 11

E-mail:

Data:

Nº:

FQ_G_13 - Rev. 01

Ref. : Upgrade de Software CIRRUS 4000 HD-OCT

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Caro Cliente Carl Zeiss,

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CARL ZEISS MEDITEC

March 25, 2011

Food and Drug Administration
San Francisco District Office
1431 Harbor Bay Parkway
Alameda, CA 94502

Dear Sir or Madam:

Enclosed please find information about a Field Corrective Action initiated by Carl Zeiss Meditec, Inc. of Dublin, California. This Field Corrective Action is for Cirrus HD-OCT Model 400 and Model 4000 with version 5.0.0 or 5.1.0 software and associated software kits.

Following the receipt of a customer complaint it was discovered that a software defect exists in versions 5.0.0 and 5.1.0 of Cirrus™ HD-OCT Instrument and Review Software. When following a specific workflow in which operators (1) acquire data, (2) log out without analyzing the data, and (3) log back in and acquire data for another patient, the exam data from the second patient may incorrectly be saved under the name of the first patient.

If you have any questions regarding this Field Corrective Action, please do not hesitate to contact me at (925) 557-4616 or at j.brimacombe@meditec.zeiss.com.

Sincerely yours,

Judith A. Brimacombe, MA
Director, Clinical/Regulatory Affairs

Enclosures: Attachment 1: Health Hazard Evaluation
 Attachment 2: Customer Notification Letter 1 with Acknowledgement Form 1
 Attachment 3: Customer Notification Letter 2 with Acknowledgement Form 2
 Attachment 4: Report of Correction 2918630-3/25/11-001-C
 Attachment 5: CD Disk containing Consignee List