Of. 16_209

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Prezado Cliente,

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
- Expert System
- Interpretation Rule Set
- + Rules
- + Patient Fields
- + Specimen Fields
- + Isolate Fields
- + Tests
- + Input
- + Export
- + OC

		Rules					
Number	Status	Automatic	Type	Description	Chart Description	Chartable	
1437	Enabled	<input checked="" type="checkbox"/>	CLSI	Isolate is resistant to one or more	Results suggest the presence of an ESBL	<input type="checkbox"/>	
1441	Enabled	<input checked="" type="checkbox"/>	CLSI	Isolate is resistant to one or more	Results suggest the presence of an ESBL	<input type="checkbox"/>	
1446	Enabled	<input checked="" type="checkbox"/>	CLSI	Isolate is a multiple-drug resistant	Isolate is a multiple-drug resistant	<input type="checkbox"/>	
1448	Enabled	<input checked="" type="checkbox"/>	CLSI	Staphylococcus aureus with a resistant	Staphylococcus aureus with a resistant	<input type="checkbox"/>	
1450	Enabled	<input checked="" type="checkbox"/>	CLSI	Oxacillin resistance (mecA) for this isolate is	Oxacillin resistance (mecA) for this isolate is	<input type="checkbox"/>	
1451	Enabled	<input checked="" type="checkbox"/>	CLSI	This isolate exhibits an antibiogram	This isolate exhibits an antibiogram	<input type="checkbox"/>	
1452	Enabled	<input checked="" type="checkbox"/>	CLSI	Isolate is a Staphylococcus spp. with	Isolate is a Staphylococcus spp. with	<input type="checkbox"/>	
1453	Enabled	<input checked="" type="checkbox"/>	CLSI	Isolate is a multiple-drug resistant	Isolate is a multiple-drug resistant	<input type="checkbox"/>	
1454	Disabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	
1456	Disabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	
1457	Disabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	
1458	Disabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	
1459	Disabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	
1460	Disabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	
1461	Disabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	
1462	Disabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	
1463	Disabled	<input checked="" type="checkbox"/>	CLSI	Isolate tested carbapenem resistant.	Isolate tested carbapenem resistant.	<input type="checkbox"/>	
1464	Enabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	
1465	Enabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	
1466	Enabled	<input checked="" type="checkbox"/>	CLSI	Isolate produces ESBL in other enhanced	Isolate produces ESBL in other enhanced	<input type="checkbox"/>	
1467	Enabled	<input checked="" type="checkbox"/>	CLSI	Isolate exhibits elevated carbapenem MICs	This isolate may be resistant to therapy with	<input type="checkbox"/>	
1468	Enabled	<input checked="" type="checkbox"/>	CLSI	The ESBL test is positive for this isolate.	The ESBL test is positive for this isolate.	<input type="checkbox"/>	
1469	Disabled	<input checked="" type="checkbox"/>	CLSI	Isolates that are resistant to all drugs tested	Isolates that are resistant to all drugs tested	<input type="checkbox"/>	
1470	Enabled	<input checked="" type="checkbox"/>	CLSI	Inducible clindamycin resistance may be	Inducible clindamycin resistance may be	<input type="checkbox"/>	
1471	Enabled	<input checked="" type="checkbox"/>	CLSI	When Klebsiella oxytoca is resistant to	When Klebsiella oxytoca is resistant to	<input type="checkbox"/>	
1472	Enabled	<input checked="" type="checkbox"/>	CLSI	When CT hyperproduction phenotype is	When CT hyperproduction phenotype is	<input type="checkbox"/>	
1473	Disabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, the SIR results for	When ESBL is detected, the SIR results for	<input type="checkbox"/>	
1474	Disabled	<input checked="" type="checkbox"/>	CLSI	When ESBL is detected, an interpretation	When ESBL is detected, an interpretation	<input type="checkbox"/>	

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BDXpert BDXpert Rules Enabled

Allow Rapid Completion of Phoenix Panels

Allow Creation of Inferred Drug Results

Rules						
Number	Status	Automatic	Type	Description	Chart Description	Chartable
101	Enabled	<input checked="" type="checkbox"/>	CLSI	Only results for ampicillin, a fluoroquinolone, and trimethoprim-sulfamethoxazole should be reported for fecal isolates of <i>Salmonella</i> spp. or <i>Shigella</i> spp. No other drugs should be reported.	Only results for ampicillin, a fluoroquinolone, and trimethoprim-sulfamethoxazole should be reported for fecal isolates of <i>Salmonella</i> spp. or <i>Shigella</i> spp. No other drugs should be reported.	<input type="checkbox"/>
102	Enabled	<input checked="" type="checkbox"/>	CLSI	Results for chloramphenicol and a third generation cephalosporin should be reported for <i>Salmonella</i> spp. from extraintestinal sites.	Results for chloramphenicol and a third generation cephalosporin should be reported for <i>Salmonella</i> spp. from extraintestinal sites.	<input checked="" type="checkbox"/>
103	Enabled	<input checked="" type="checkbox"/>	CLSI	Results for chloramphenicol and a third generation cephalosporin should be reported for <i>Salmonella</i> spp. from extraintestinal sites.	Results for chloramphenicol and a third generation cephalosporin should be reported for <i>Salmonella</i> spp. from extraintestinal sites.	<input type="checkbox"/>
104	Enabled	<input checked="" type="checkbox"/>	CLSI	<i>Salmonella</i> spp. and <i>Shigella</i> spp. should not be tested for first and second generation cephalosporins even if they appear to be susceptible in vitro because they are not effective clinically.	<i>Salmonella</i> spp. and <i>Shigella</i> spp. should not be tested for first and second generation cephalosporins because they are not effective clinically.	<input checked="" type="checkbox"/>
105	Enabled	<input type="checkbox"/>	CLSI	Ampicillin is the class.	Ampicillin is the class.	<input type="checkbox"/>

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