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Urgent Field Safety Notice Attachment

This notice must be passed on to all those who need to be aware within your organisation or to any organisation where the devices have been transferred to.

**Eckert & Ziegler
BEBIG GmbH**


Robert-Rössle-Str. 10
13125 Berlin
Germany
www.bebig.com

Subject: Offset of Depth Dose Curve
Reference: **UFSN-2019Ru106**
Date of notification: 15 Aug 2019

With this attachment, we provide to you the corrected dose rate curve $\dot{D}(z)$ along the central axis of each individual Ru-106 Eye Applicator concerned by the offset. The correction has been performed by interpolation to the dose rate at the respective distances to the surface z_{new} .

The graph replaces that one given on the appendix of the certificate. In the table, the measured dose rates \dot{D}_{meas} at the corrected distances ($z_{\text{corr}} = 1.35, 2.35, \dots, 10.35$ mm) are stated together with the interpolated data \dot{D}_{int} ($z_{\text{new}} = 1.00, 2.00, \dots, 10.00$ mm).

Serial number of applicator CCB 2735
Certificate no. 51119
Reference date 26.04.2019

* Due to technical reasons the value at 0.83 mm was measured correctly and therefore does not change.
** Extrapolated data according the previous and new fit function.
*** Relative deviation between previous and corrected dose rates

Please note: Since we have not carried out a new measurement, we are unfortunately unable to issue a new certificate for you. Hence please use this information as a substitute for the corresponding appendix of the certificate.

Managing Directors:
Dr. Harald Hasselmann
Dirk Warmuth

Berlin Office:
Amtsgericht Charlottenburg
Reg.-Nr. HRB 42949 B
USt.-ID-Nr.: DE 137169788

Bank Account:
Commerzbank AG
Kto.: 0424648
BLZ: 12040000

Urgent Field Safety Notice Attachment

This notice must be passed on to all those who need to be aware within your organisation or to any organisation where the devices have been transferred to.

**Eckert & Ziegler
BEBIG GmbH**


Robert-Rössle-Str. 10
13125 Berlin
Germany
www.bebig.com

Subject: Offset of Depth Dose Curve
Reference: **UFSN-2019Ru106**
Date of notification: 15 Aug 2019

With this attachment, we provide to you the corrected dose rate curve $\dot{D}(z)$ along the central axis of each individual Ru-106 Eye Applicator concerned by the offset. The correction has been performed by interpolation to the dose rate at the respective distances to the surface z_{new} .

The graph replaces that one given on the appendix of the certificate. In the table, the measured dose rates \dot{D}_{meas} at the corrected distances ($z_{corr} = 1.35, 2.35, \dots, 10.35$ mm) are stated together with the interpolated data \dot{D}_{int} ($z_{new} = 1.00, 2.00, \dots, 10.00$ mm).

Serial number of applicator: COB 1252
Certificate no.: 51219
Reference date: 26.04.2019

* Due to technical reasons the value at 0.83 mm was measured correctly and therefore does not change.
** Extrapolated data according the previous and new fit function.
*** Relative deviation between previous and corrected dose rates

Please note: Since we have not carried out a new measurement, we are unfortunately unable to issue a new certificate for you. Hence please use this information as a substitute for the corresponding appendix of the certificate.

Managing Directors:
Dr. Harald Hasselmann
Dirk Warmuth

Berlin Office:
Amtsgericht Charlottenburg
Reg.-Nr. HRB 42949 B
USt.-ID-Nr.: DE 137169788

Bank Account:
Commerzbank AG
Kto.: 0424648
BLZ: 12040000