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⚠ **Warning:** Do not send a patient home if:

- there is a pump memory error.
- there is an unresolved alarm condition.
- telemetry was interrupted during pump update or interrogation.
- pump status was not confirmed.

Failure to correct pump memory errors, alarm conditions, or incomplete telemetry can result in incorrect or incomplete programming, unknown or unexpected pump status, and a clinically significant or fatal drug underdose or overdose.

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7. When the programming head is extended, a status light is visible on the back of the programming head:

- **Green light flashing**—Telemetry is in progress
- **Amber light flashing**—Telemetry has not been established

⚠ **Caution:** To ensure successful telemetry, hold the programming head steady over the implanted pump until confirmation is received that telemetry was successful. If the programming head is not held over the pump properly, telemetry may not be established or may be interrupted, which can result in incorrect or incomplete programming.

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#### Telemetry failure checklist

The most common corrections for telemetry failures are listed below. Refer to Table 6 on page 77 for explanations of specific error messages.

- Ensure the programming head is undocked and the cable is fully extended.
- Decrease the distance between the programming head and the pump by pressing the programming head firmly over the implanted pump.
- Hold the programming head steady over the pump, then press the PROGRAMMING key or select the PROGRAM button.
- If telemetry fails while holding the programming head steady, move the programming head around slowly near the implanted pump. Press the PROGRAMMING key, or select the PROGRAM button.
- When holding the programming head over the implanted pump, do not drape your hand over the back of the programming head. Hold the programming head at the base.
- Ensure that the programmer batteries are not low or depleted.
- Move away from sources of possible EMI (eg, computer monitor, cell phone, motorized wheelchair).

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
- ⚠ Warning:** Do not send a patient home if:
- there is a pump memory error.
  - there is an unresolved alarm condition.
  - telemetry was interrupted during pump update or interrogation.
  - pump status was not confirmed.

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**⚠ Warning:** To ensure successful telemetry, hold the communicator steady over the implanted pump until confirmation is received that telemetry was successful. If the communicator is not held over the pump properly, telemetry may not be established or may be interrupted, which can result in unintended or incomplete programming.

#### **If telemetry is interrupted during a pump update**

When telemetry is interrupted, either by user action or system action, the pump state is unknown, and the pump may or may not be stopped. The system responds to interrupted telemetry with a high-level alert.

To continue updating the pump, tap the **CONTINUE UPDATE** button. If the application is closed and the pump is left in an unknown state, start a programming session with the pump and check all of the current pump settings.

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