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Matriz Lagoa Santa  
Rua Alfredo Albano da Costa, 100  
Sls. 1, 2 e 3 - Dist. Ind. Genesco A. de Oliveira  
Lagoa Santa - MG  
Brasil - 33.400-000  
Tel.: +55 31 3689 6600

Filial São Paulo  
Av. Dr. Chucri Zaidan, 1.240  
19º Andar - Conj. 1.902 e 1.904 - Torre Golden  
Ed. Morumbi Corporate - Vila São Francisco  
São Paulo - SP  
Brasil - 04.709-111  
Tel.: +55 11 3065 7550

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Caro cliente

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### Subtela Setup/Test/Data Setup

The screenshot displays the 'Data Setup' configuration interface. At the top, it shows the current test as 'NeoHb', the reagent set as '90085-A', and the method version as '1.016'. There are 'Update Kit' and 'Edit Test' buttons. Below this, there are two main sections: 'Reagents / Cartridges Log' with radio buttons for 'Cartridges' and 'Reagents', and 'Settings' with radio buttons for 'Peak Table', 'Integration Table', 'Data Setup' (which is selected), and 'Pattern Setting'. An 'Allow Editing' button is located below the settings. The main area is titled 'User Selected Sample Types' and contains a list of checkboxes: 'Blank' (checked), 'Primer' (checked), 'FAES control' (checked), 'FADC control' (checked), 'Tray control' (unchecked), 'Proficiency control' (unchecked), 'Positional control' (unchecked), and 'Carryover control' (unchecked). At the top right, there is a 'Select New Test' dropdown menu set to 'NeoHb'.

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## Subtela Setup/Test/Pattern Setting

The screenshot displays the 'Pattern Setting' configuration window. At the top, it shows 'Current Test: NeoHb', 'Reagent Set: 90685-A', and 'Method Version: 1.016'. There are buttons for 'Update Kit' and 'Edit Test'. Below this, there are two main sections: 'Reagents / Cartridges Log' with radio buttons for 'Cartridges' and 'Reagents', and 'Settings' with radio buttons for 'Peak Table', 'Integration Table', 'Data Setup', and 'Pattern Setting' (which is selected). An 'Allow Editing' button is located below the settings. The main configuration area includes several input fields and checkboxes:

- FAST Area Threshold (%):** From: 24.0, To: 24.0
- Total Area Limit (AU):** From: 750000, To: 6300000
- FAST Height / Total Area threshold:** 0.0300
- F peak Height / Area Threshold for FADC:** 0.22
- F Only Total Area Lower Limit:** 150000
- Minimum Area% of Peak Cutoff A:** 0
- Checkboxes:**
  - Display pattern in data views
  - Display pattern in reports
  - Add an asterisk at the start of the pattern
  - If the A peak is less than or equal to its Min Area% threshold
  - If a peak is omitted from the pattern due to the One Fourth Rule
  - If the FAST Height-to-total-area ratio is above its threshold
  - Enable the One Fourth Rule

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