## **URGENTE: RECALL DE DISPOSITIVO MÉDICO**

Cassetes VERITAS™ Advanced Infusion (VRT-AI) e Cassetes VERITAS™ Advanced Fluidics (VRT-AF)

01 de Junho de 2023

Prezado cliente da Johnson & Johnson Surgical Vision:

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Data: \_\_\_\_\_

Exemplo de etiqueta da tampa do Cassete VRT-AI

The diagram shows a purple-bordered label for the VRT-AI Advanced Infusion Pack. On the left, a vertical purple bar contains the text 'VRT-AI' and 'Advanced Infusion Pack'. The main label area features the 'veritas' logo in a large, grey, lowercase font, with 'Advanced Infusion' in a smaller font below it. To the left of the main text, it lists the contents: 'Contains (1) each: 1. Cassette and Tubing Assembly, 2. Monitor Drape Cover, 3. Mayo Stand Drape Cover, 4. Test Chamber'. The Johnson & Johnson VISION logo is centered, with the company address: 'Johnson & Johnson Surgical Vision, Inc., 1700 E. St. Andrew Place, Santa Ana, CA 92705 USA'. Below this, it states 'Product of Mexico'. To the right, there are regulatory symbols: 'Rx only', 'STERILE EO', 'MD', and 'CE 0413'. A 'REF' box contains 'VRT-AI', and an 'EC REP' box contains 'AMO Ireland, Block B, Liffey Valley Office Campus, Quarryvale, Co. Dublin, Ireland'. At the bottom, there are two barcode areas. The first is a standard barcode with the number '(01)050504 74700901'. The second is a larger barcode with the number '(17)YYMMDD(10)TEST1234'. Above the second barcode, there are two date fields: 'YYYY-MM-DD' and 'YYYY-MM-DD'. A red box highlights the 'REF' and 'EC REP' information, with a callout pointing to the 'REF' text. Another red box highlights the 'LOT' and 'TEST1234' information, with a callout pointing to the 'TEST1234' text. The label also includes 'PATENTS: www.jnjvisionpro.com/patents', '© Johnson & Johnson Surgical Vision, Inc. 2021', and 'Z353726 Rev. E 02/2021'.

VRT-AI  
Advanced Infusion Pack

veritas  
Advanced Infusion

Contains (1) each:  
1. Cassette and Tubing Assembly  
2. Monitor Drape Cover  
3. Mayo Stand Drape Cover  
4. Test Chamber

Johnson & Johnson VISION  
Johnson & Johnson Surgical Vision, Inc.  
1700 E. St. Andrew Place  
Santa Ana, CA 92705 USA  
Product of Mexico

Rx only  
STERILE EO  
MD  
www.e-ft.com  
CE 0413

REF VRT-AI  
EC REP AMO Ireland  
Block B  
Liffey Valley Office Campus  
Quarryvale, Co. Dublin, Ireland

PATENTS: www.jnjvisionpro.com/patents  
© Johnson & Johnson Surgical Vision, Inc. 2021  
Z353726 Rev. E 02/2021

LOT TEST1234  
YYYY-MM-DD  
YYYY-MM-DD  
(01)050504 74700901  
(17)YYMMDD(10)TEST1234 Z353726 Rev. E 0520

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Exemplo:  
Local do  
Número do  
Lote

## Exemplo de etiqueta da tampa do Cassete VRT-AF

The diagram shows a rectangular label for the VRT-AF cassette cover. On the left side, there is a vertical blue bar with the text "VRT-AF" and "Advanced Fluidics Pack". The main label area contains the following information:

- veritas** (large font) and **Advanced Fluidics** (smaller font)
- Contains (1) each:
  - 1. Cassette and Tubing Assembly
  - 2. Monitor Drape Cover
  - 3. Mayo Stand Drape Cover
  - 4. Test Chamber
- Johnson & Johnson VISION** logo and address: Johnson & Johnson Surgical Vision, Inc., 1700 E. St. Andrew Place, Santa Ana, CA 92705 USA
- Product of Mexico
- REF** VRT-AF (highlighted with a red box and callout)
- EC REP** AMO Ireland, Block B, Liffey Valley Office Campus, Quarryvale, Co. Dublin, Ireland
- STERILE EO** and **MD** symbols
- CE** 0413 symbol
- LOT** TEST 1234 (highlighted with a red box and callout)
- Barcode with lot number TEST 1234
- Patents: [www.jnjvisionpro.com/patents](http://www.jnjvisionpro.com/patents)
- © Johnson & Johnson Surgical Vision, Inc. 2021
- Z353725 Rev. D 02/2021

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Exemplo:  
Local do  
Número do  
Lote